IVF Basics® SpermiTec M

EN

USED ABBREVIATIONS

IUI IntraUterine Insemination
HSA Human Serum Albumin
LAL Limulus Amebocyte Lysate
SSCP Summary of Safety and Clinical Performance
MSDS Material Safety Data Sheet

GENERAL

IVF Basics® SpermiTec M has been developed for all human sperm washing and sperm preparation techniques.

IVF Basics® SpermiTec M contains HEPES, no CO₂ incubation is required.

APPLICATION

Intended use

IVF Basics® SpermiTec M is used for all human sperm washing and sperm preparation techniques including swim-up and density gradient centrifugation. It is also used for Intrauterine insemination (IUI).

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

IVF Basics® SpermiTec M is a ready-to-use HEPES-buffered medium which contains bicarbonate, physiologic salts, glucose, lactate, pyruvate, phenol red, 10mg/liter gentamicin sulphate (medicinal substance) and Splitter HSA (medicinal substance derived from human blood plasma).

The MSDS for IVF Basics® SpermiTec M also describes the product composition and is available on the website of Gynotec B.V.

QUALITY CONTROL

- pH: 7.20-7.60
- Osmolality: 270-290 mOsm/kg
- Endotoxins test by LAL methodology (USP <85>): < 0.25 EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96h): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic techniques should be used to avoid possible contamination, even when the products contain gentamicin.
- Always wear protective clothing when handling specimens.
- Handle specimens as if capable of transmitting HIV or hepatitis.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of human blood products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured in European Pharmacopoeia specifications by established processes.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

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