APPLICATION

SpermTec® G-80 and G-45 gradients are ready to use gradient systems for semen preparation. The SpermTec® gradient system can be used for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF) and IntraCyttoplasmatic Sperm Injection (ICSI).

For professional use only.

COMPOSITION

SpermTec® gradients consist of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution) supplemented with phenol red, 0.4-2.2g/Albumina Humana (albumin (medical substance derived from human blood plasma) and 10mg/l gentamicin (medicinal substance).

QUALITY CONTROL

- pH between 7.20-7.90 (Release criteria: 7.20-7.60).
- Osmolality:
 - G-80: 290-330 mOsm/kg
 - G-45: 280-310 mOsm/kg
- Density:
 - G-80: 1.097-1.107 g/ml
 - G-45: 1.0500-1.0700 g/ml
- Endotoxin (USP <85>): < 0.5EU/ml
- Sterility test according to current Ph. Eur. 2.6.1/ USP <71>. No growth
- Human sperm survival assay:
 - % motility compared with control after 4 hours: $\geq 80\%$
 - % motility compared with control after 24 hours: $\geq 75\%$
- Chemical composition
- Use of P1, USP grade products if applicable
- Not MEA tested
- Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

PRECAUTIONS AND WARNINGS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for product exclusion. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Gynotec B.V. and, if applicable, to the authority of the EU Member State in which the user and/or patient is established.

PRE-ISE CHECKS

- Do not use the product if the seal of the container is opened or defect when the product is delivered.
- Do no use if the products shows any evidence of microbial contamination or becomes cloudy.
- Do not use after expiry date.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium and avoid exposure to light.
- Gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.

STORAGE CONDITIONS

- Store products between 2-8 °C.
- Keep away from (sun)light.
- After opening the container, do not use the product more than 7 days. Sterile conditions must be maintained and product must be stored at 2°C – 8°C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

$g = 1.118 \times r \times rpm^2$		$rpm = \sqrt{\frac{g}{1.118 \times r}}$
-----------------------------------	---	---

<i>r = radius of centrifuge in mm</i> <i>rpm = rotation per minute / 1000</i>		
--	--	--

Example 1 <i>r</i> = 100 mm <i>rpm</i> = 1800 rotations per minute <i>g</i> = 1.118 x 100 x 324 = 362 g	Example 2 <i>r</i> = 100 mm <i>g</i> = 350g <i>rpm</i> = $\sqrt{350 / (1.118 \times 100)}$ = 1.77 <i>rpm</i> = 1770 rotations per minute
---	---

INSTRUCTIONS FOR USE

Method
Each laboratory should establish its own validated procedures and protocols.

- Fresh semen samples**
- Mix the density gradient bottles by 5 bottle inversions before use.
 - Example 1**
r = 100 mm
rpm = 1800 rotations per minute
g = 1.118 x 100 x 324 = 362 g
 - Example 2**
r = 100 mm
g = 350g
rpm = $\sqrt{350 / (1.118 \times 100)}$ = 1.77
rpm = 1770 rotations per minute
 - Bring all components of the system and samples to room temperature or to 37 °C.
 - Transfer 2.5ml of SpermTec® G-45 into a sterile disposable centrifuge tube.
 - Place 2.5ml of SpermTec® G-80 under the SpermTec® G-45 layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the SpermTec® G-80. This two layer gradient is stable for up to two hours.
 - Gently place up to 2.5ml of liquefied semen onto the SpermTec® G-45 using a transfer pipette. Do not use a high volume than the volumes of the individual gradient layers or more than 10⁶ cells.
 - Centrifuge for 15 to 18 minutes at 350g to 400g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
 - Remove supernatant down to the pellet.
 - Add 2-3ml of sperm washing medium (e.g. SpermTec® Wash or SpermWash® distributed by Gynotec B.V.) and resuspend the pellet.
 - Centrifuge for 8 to 10 minutes at 300g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.
 - Remove supernatant down to the pellet and repeat steps 7 and 8.
 - Remove supernatant and replace with a suitable volume of appropriate medium.

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SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)
The SSCP for SpermTec® G-80 and G-45 describes safety and performance characteristics for the method of sperm selection stable on the website of Gynotec B.V. (www.gynotec.nl).

For further questions regarding to the safety and performance, please contact Gynotec B.V. for customer or technical support.

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