**IVF Basics® HTF**  
**IVF Basics® HTF 0.4% HSA**  
-

**IVF Basics® HTF HEPES**  
**IVF Basics® HTF HEPES 0.4% HSA**

**IVF Basics® HTF HEPES**  
**IVF Basics® HTF HEPES 0.4% HSA**

## EN

### USED ABBREVIATIONS

ICSI Intracytoplasmatic Sperm Injection  
IVF In Vitro Fertilization  
IUI IntraUterine Insemination  
HSA Human Serum Albumin  
LAL Limulus Amebocyte Lysate  
SSCP Summary of Safety and Clinical Performance  
MSDS Material Safety Data Sheet

### GENERAL

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) have been developed for washing, handling, conditioning and culturing of human gametes and embryos.

IVF Basics® HTF (0.4% HSA) should be pre-equilibrated overnight or at least for a minimum of 4 hours before use at 37°C and 5% CO<sub>2</sub>.

IVF Basics® HTF HEPES (0.4% HSA) contains HEPES, no CO<sub>2</sub> incubation is required.

### APPLICATION

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) are designed for washing and handling of human gametes, and for handling and conditioning of human embryos. Also is used for embryo transfer.

### For professional use only.

IVF Basics® HTF (0.4% HSA) can be used for the following procedures:

- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling of human embryos
- IUI, IVF, ICSI
- Embryo culture from day 1 to blastocyst stage
- Embryo transfer

IVF Basics® HTF HEPES (0.4% HSA) can be used for the following procedures:

- Flushing for oocyte collection
- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling of human embryos
- IUI, ICSI
- Embryo transfer

**COMPOSITION**  
IVF Basics® HTF and IVF Basics® HTF HEPES consist of a balanced physiologic salt solution, supplemented with carbohydrate energy sources such as glucose, pyruvate and lactate, and 10mM gentamicin sulphate (medicinal substance). IVF Basics® HTF is bicarbonate buffered, while IVF Basics® HTF HEPES additionally contains HEPES.

It is strongly suggested to supplement IVF Basics® HTF and IVF Basics® HTF HEPES with 4g/liter HSA (medicinal substance derived from human blood plasma) before use.

IVF Basics® HTF 0.4% HSA and IVF Basics® HTF HEPES 0.4% HSA are presupplemented with 4g/liter HSA (medicinal substance derived from human blood plasma) and are ready-to-use.

The MSDS for IVF Basics® HTF media also describes the product composition and is available on the website of Gynotec B.V.

**MATERIAL NOT INCLUDED**  
IVF Basics® HTF and IVF Basics® HTF HEPES (media without HSA): pharmaceutical grade Human Albumin Solution with registered Plasma Master File

**QUALITY CONTROL**  
pH: 7.20-7.60 (± 0.5% CO<sub>2</sub> for IVF Basics® HTF (with HSA))  
Osmolarity: 270-290 mOsm/kg

Endotoxin test by LAL methodology (USP <85>: < 0.25 EU/ml)  
Sterility test by the current Eur. Ph. 2.6.1/1 USP <71>: No growth

One-cell mouse embryo assay (blastocysts after 96h): ≥ 80%  
Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80%

Chemical composition

Use of Ph Eur or USP grade products if applicable  
Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

### PRECAUTIONS AND WARNINGS

Aseptic techniques should be used to avoid possible contamination, even when the products contain gentamicin.

Always wear protective clothing when working with specimens.

Handle specimens as if capable of transmitting HIV or hepatitis.

All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

Standard measures to prevent infections resulting from the use of medical products prepared from human blood or plasma must include selection of appropriate individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products are prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

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Always wear protective clothing when working with specimens.

Handle specimens as if capable of transmitting HIV or hepatitis.

All blood products should be treated as potentially infectious.

Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. Known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

Standard measures to prevent infections resulting from the use of medical products prepared from human blood or plasma must include selection of appropriate individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses.

Despite this, when medicinal products are prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens.

There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

For further questions regarding to the safety and performance, please contact Gynotec B.V.

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