IVF Basics® HTF
IVF Basics® HTF 0.4% HSA

IVF Basics® HTF HEPES
IVF Basics® HTF HEPES 0.4% HSA

**EN**

**USCI ABBREVIATIONS**

ICSI Intracytoplasmatische Sperm Injection
IVF In Vitro Fertilization
IUI Intra uterine insemination
HSA Human Serum Albumin
LAL Limulus Amebocyte Lystate
SSCP Summary of Safety and Clinical Performance
MSDS Material Safety Data Sheet

**GENERAL**

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) have been developed for washing, handling, conditioning and culturing of human gametes and embryos.

- IVF Basics® HTF (0.4% HSA) should be pre-equilibrated overnight or at least for a minimum of 4 hours before use at 37°C and 5% CO<sub>2</sub>.
- IVF Basics® HTF HEPES (0.4% HSA) contains HEPES, no CO<sub>2</sub> incubation is required.

**APPLICATION**

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) are designed for washing and handling of human gametes, and for handling and conditioning of human embryos. Also to be used for embryo transfer.

For professional use only.

IVF Basics® HTF (0.4% HSA) can be used for the following procedures:

- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, IVF, ICSI
- Embryo culture from day 1 to blastocyst stage
- Embryo transfer

IVF Basics® HTF (0.4% HSA) can be used for the following procedures:

- Flushing for oocyte collection
- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, ICSI
- Embryo transfer

IVF Basics® HTF HEPES (0.4% HSA) can be used for the following procedures:

- Flushing for oocyte collection
- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, ICSI
- Embryo transfer

IVF Basics® HTF HEPES (0.4% HSA) can be used for the following procedures:

- Flushing for oocyte collection
- Washing/handling of human ova
- Washing/handling of spermatozoa
- Swim-up of spermatozoa
- Production of density gradient media
- Washing/handling for human embryos
- IUI, ICSI
- Embryo transfer

IVF Basics® HTF (0.4% HSA) and IVF Basics® HTF HEPES (0.4% HSA) consist of a balanced physiologic salt solution, supplemented with carbohydrate energy sources such as glucose, pyruvate and lactate, and 10mg/liter gentamicin sulphate (medicinal substance). IVF Basics® HTF is bicarbonate buffered, while IVF Basics® HTF HEPES additionally contains HEPES.

It is strongly suggested to supplement IVF Basics® HTF and IVF Basics® HTF HEPES with Agfiter HSA (medicinal substance derived from human blood plasma) before use.

IVF Basics® HTF 0.4% HSA and IVF Basics® HTF HEPES 0.4% HSA are presupplemented with 4g/liter HSA (medicinal substance derived from human blood plasma) and are ready-to-use.

The MSDS for IVF Basics® HTF media also describes the product composition and is available on the website of Gynotec B.V.

**MATERIAL NOT INCLUDED**

IVF Basics® HTF and IVF Basics® HTF HEPES (media without HSA): pharmaceutical grade Human Albumin Solution with registered Plasma Master File

**QUALITY CONTROL**

- pH: 7.20-7.60 (At 5% CO<sub>2</sub> for IVF Basics® HTF (with HSA))
- Osmolality: 270-290 mOsm/kg
- Endotoxin test by LAL methodology (USP <85>): < 0.25 EU/ml
- Sterility test by the current Ph. Eur. 2.6.1 USP <71>: No growth
- Sterility test by the current Ph. Eur. 2.6.1 USP <71>: Geen groei
- Human sperm survival assay (% motility compared with control after 24 hours): ≥ 80%
- Chemical composition
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request

Sterilized by using aseptic processing techniques.

**PRECAUTIONS AND WARNINGS**

Aseptic techniques should be used to avoid possible contamination, even when the products contain gentamicin.

- Always wear protective clothing when working with specimens.
- Place the syringe onto a catheter and eject the medium.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-evaluated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

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► Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro N.V. and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

**PRE-USE CHECKS**

- Do not use the product if:
  - the seal on the bottle is broken or open when the product is delivered; or it becomes discoloured; shows any evidence of microbial contamination or cloudy;
  - expiry date has been exceeded.
- Do not freeze before use.
- Do not resterilize after opening.
- Products which include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.

► Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of medium. Discard excess (unused) media.

**STORAGE CONDITIONS**

- Store refrigerated (2°C – 8°C).
- Keep away from (sun)light.
- After opening the container, do not use the product longer than 7 days. Sterile conditions must be maintained and product must be stored at 2°C – 8°C.
- Discard the devices in accordance with local regulations for disposal of medical devices.

**INSTRUCTIONS FOR USE**

Each laboratory should establish its own validated procedures and protocols.

Addition of human albumin solution before use:
It is strongly recommended to supplement IVF Basics® HTF media with human albumin with 4 g/l HSA (see material not included). Use aseptic techniques to avoid introduction of microorganisms to the medium.

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- Plaats de syringe op de katheter en spuit het medium.
- Suck up the IVF Basics® HTF media into the syringe and push down the piston to the 10µl calibration mark.
- Suck up the embryo(s) into the catheter in such a way that the volume is approximately 20-30µl.
- Hand over catheter and syringe to the clinician for insertion and transfer into the uterus.

**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)**
The SSCP for IVF Basics® HTF media describes safety and performance characteristics for the media and is available on the website of Gynotec B.V. (www.gynotec.nl).

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Catalogusnummer

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IVF Basics® HTF
IVF Basics® HTF 0.4% HSA

Product code	Product description	Product code	Product description
HTF-B-10	1 x 10 ml of IVF Basics® HTF	HTF-B-0.4A-10	1 x 10 ml of IVF Basics® HTF 0.4% HSA
HTF-B-30	1 x 30 ml of IVF Basics® HTF	HTF-B-0.4A-30	1 x 30 ml of IVF Basics® HTF 0.4% HSA
HTF-B-60	1 x 60 ml of IVF Basics® HTF	HTF-B-0.4A-60	1 x 60 ml of IVF Basics® HTF 0.4% HSA
HTF-B-120	1 x 120 ml of IVF Basics® HTF	HTF-B-0.4A-120	1 x 120 ml of IVF Basics® HTF 0.4% HSA
HTF-B-250	1 x 250 ml of IVF Basics® HTF	HTF-B-0.4A-250	1 x 250 ml of IVF Basics® HTF 0.4% HSA
HTF-B-500	1 x 500 ml of IVF Basics® HTF	HTF-B-0.4A-500	1 x 500 ml of IVF Basics® HTF 0.4% HSA

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