



Certificate of Analysis

Product: SpermFilter® 100%

Batch number: FLT100-812
Expiration Date: 30-11-2019
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,31
Osmolality:	300 – 330 mOsm/kg	305
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
Endotoxin:	< 0,5 EU/ml	0,174 EU/ml
Density:	1,1150 – 1,1250 g/ml	1,1197
Sperm Survival:	≥ 80%	98%

Approved on:
(DD/MM/YYYY)

By: J. Berbers

Signature:

Quality Control Department, Gynotec B. V.

25-05-2018